**VALEDA®**

Light Delivery System

User Manual

# Indications for Use

The indicated use is for treatment of ocular damage and disease using photobiomodulation, including inhibition of inflammatory mediators, oedema or drusen deposition, improvement of wound healing following ocular trauma or surgery, and increase in visual acuity and contrast sensitivity in patients with degenerative diseases such as dry age-related macular degeneration.

## Contents

- 1. Introduction ..... 1
- 2. Valeda Light Delivery System Description ..... 1
- 3. Valeda Light Delivery System Operation ..... 3
- 4. Labelling Symbols .....4
- 5. Indications for Use ..... 5
- 6. Contraindications for Use ..... 5
- 7. Important Safety Instructions ..... 6
- 8. Installation ..... 7
- 9. Operating Instructions ..... 7
- 10. Valeda Treatment Credits ..... 17
- 11. Valeda Settings ..... 19
- 12. Cleaning ..... 20
- 13. Disposal ..... 20
- 14. Maintenance ..... 20
- 15. Troubleshooting ..... 21
- 16. Technical Specifications ..... 23
- 17. Glossary ..... 24
- 18. Guidance and Manufacturer’s Declarations ..... 25

## 1. Introduction

This manual provides the indications, contraindications, instructions for use, warnings, cautions and precautions for the Valeda® Light Delivery System. Carefully read this manual in its entirety before using the Valeda Light Delivery System. Failure to follow these instructions may result in improper use of the system.

Terminology – within this User Manual, the words “Valeda” and “System” are used to identify collectively all components (and accessories) of the Valeda Light Delivery System. For questions regarding Valeda, please contact LumiThera®.

## 2. Valeda Light Delivery System Description

### Overview

Valeda is a multi-wavelength light emitting diode (LED) system designed for eye care professionals to use in the treatment of the eye with photobiomodulation (PBM). Valeda will provide a pre-set treatment of PBM to the patient’s eye and retinal tissue through the open and closed eyelid.

Valeda uses three LEDs to generate 590, 660 and 850 nm light, which is conditioned through a series of optics to produce a uniform, noncoherent beam of 30 mm diameter at the plane of treatment. The beam is directed by adjustable mirrors to allow treatment of either eye without repositioning the patient.

The operator interface (Figure 1) consists of a touchscreen, start and stop push buttons, and a joystick to centre the beam on the patient’s eye. A USB port is used to load Valeda Treatment Credits into the system.

The patient interface (Figure 2) consists of a fixed forehead rest, the light aperture and an adjustable chin rest. The adjustment range of the chin rest is designed to accommodate the full range of the intended patient population. Single-use disposable protective barriers mitigate contamination between patients.

FIGURE 1 OPERATOR'S PERSPECTIVE

FIGURE 2 PATIENT'S PERSPECTIVE

**Components**

Description Reference	Reference
Valeda Light Delivery System	200001
Valeda User Manual	LBL-0001
Valeda Chin Rest Protective Barriers	300001
Valeda Forehead Rest Protective Barriers	300002
Valeda LumiKey	300003
Valeda Dust Cover	300004
Valeda Power Cord(s)	300005
Valeda Treatment Credits	300006

Replacement components may be purchased by contacting your LumiThera representative or through the LumiThera e-commerce store at [www.lumithera.com](http://www.lumithera.com).

## Disposables

There are 2 types of disposable protective barriers: one for the forehead rest and one for the chin rest (Figure 3).

### Valeda LumiKey

A Valeda LumiKey is a standard USB storage drive used to transfer Valeda Treatment Credits as electronic files (LumiFiles) generated by LumiThera to Valeda. Encrypted electronic files ensure that Valeda will process only files prepared by LumiThera.

### Valeda Dust Cover

A dust cover is provided to protect Valeda from dust and mild abrasions. Remove before and replace after use of Valeda.

## 3. Valeda Light Delivery System Operation

1. The patient is seated in front of Valeda and disposable protective barriers are placed over the forehead and chin rests between Valeda and the patient.
2. The patient rests against Valeda forehead rest for stabilisation.
3. The operator adjusts the chin rest to best maintain the patient's eye position relative to Valeda, prompts the patient to open or close their eyes as necessary, performs final alignment of the light source to the patient's eye, activates the light source, and monitors the treatment progress.
4. During treatment, the patient positions their head against Valeda, maintains that position during treatment (35-90 seconds for each delivery interval), and opens or closes their eyelids as instructed by the operator.
5. The patient does not have to fixate their eye position or focus during light delivery.
6. Valeda will indicate when the LEDs are on and automatically turn off the light source after the defined treatment time.
7. The patient then pulls away from Valeda, the protective barriers are discarded and Valeda is turned off.
8. The system does not require cleaning between treatments, but may be cleaned as described in this manual. The dust cover may be placed over the system when appropriate.

## 4. Labelling Symbols

	<b>Consult instructions for use.</b> Indicates the need for the user to consult the instructions for use.
	<b>Refer to instruction manual/booklet.</b> Signifies that the instruction manual/booklet must be read.
	<b>Manufacturer.</b> Indicates the medical device manufacturer.
	<b>Date of manufacture.</b> Indicates the date when the medical device was manufactured.
	<b>Authorised Representative.</b> Indicates the authorised representative in the European Community.
	<b>Catalogue number.</b> Indicates the manufacturer's catalogue number so that the device can be identified.
	<b>Batch code.</b> Indicates the manufacturer's batch code so that the batch or lot can be identified.
	<b>Serial number.</b> Indicates the manufacturer's serial number so that the device can be identified.
	<b>Applied Part, Type BF.</b>
	<b>Do not reuse.</b> Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.
	<b>Temperature limit.</b> Indicates the temperature limits to which the medical device can be safely exposed.
	<b>Humidity limitation.</b> Indicates the range of humidity to which the medical device can be safely exposed.
	<b>WEEE Directive.</b> Waste Electrical and Electronic Equipment Directive.
	<b>Fragile, handle with care.</b> Indicates a medical device that can be broken or damaged if not handled carefully.
	<b>Keep dry.</b>
	<b>Do not stack.</b>
	<b>This way up.</b>

	<b>Recycle.</b>
	<b>Atmospheric pressure limitation.</b> Indicates the range of atmospheric pressure to which the medical device can be safely exposed.
	<b>Packaging unit.</b> Indicates the number of pieces in the package.
	<b>Alternating current.</b> Indicates on the rating plate that the equipment is suitable for alternating current only.
	<b>Fuse.</b> Identifies fuse boxes or their location.
	<b>Valeda Treatment Credits.</b>
	<b>Back.</b> Display navigation button to return to the prior screen.
	<b>Continue.</b> Display navigation button to proceed to the next screen.
	<b>Confirmation.</b> Display button to confirm displayed information.
	<b>Settings.</b> Display navigation button to proceed to the system settings.
	<b>Home.</b> Display navigation button to return to the home screen.
	<b>Information.</b> Display navigation button to proceed to the system information screen.

## 5. Indications for Use

The Valeda Light Delivery System is indicated for the treatment of ocular damage and disease using photobiomodulation, including inhibition of inflammatory mediators, oedema or drusen deposition, improvement of wound healing following ocular trauma or surgery, and increase in visual acuity and contrast sensitivity in patients with degenerative diseases such as dry age-related macular degeneration.

## 6. Contraindications for Use

As a precaution, patients have not been tested and should not be treated with Valeda if they have any known photosensitivity to yellow light, red light or near infrared radiation (NIR), or if they have a history of light activated central nervous system disorders (e.g., epilepsy, migraine). In addition, patients should not receive treatment within 30 days of using any photosensitising agent (e.g., topicals, injectables) before consulting with their physician.

## 7. Important Safety Instructions

Valeda is a medical device and should be operated only by qualified medical professionals in accordance with this manual and any applicable country or local laws and regulations. Prior to using Valeda, be sure to review all described warnings and entire instructions for use. Failure to follow the warnings may result in injury to the patient or operator or cause property damage.

It is the responsibility of the operator and the institution or organisation where Valeda will be operated to comply with their safety standards concerning the use of this system.

### Warnings

Warnings are information relating to potential dangers when using Valeda.

- Do not use the system if damage is observed or if the system is not operating as expected.
- Do not modify the system.
- Use only the power cord provided with the system or supplied by LumiThera.
- Do not use the power cord if it has been damaged. Use of a damaged power cord may result in a fire or electrical shock hazard.
- To avoid the risk of electrical shock, the system must be connected only to supply mains with protective earth.
- Do not treat patients with open sores that may come into contact with the system.
- Do not use this system on any part of the body other than designated treatment sites.
- Disconnect power to the system before cleaning.
- Patients should remove their glasses and contact lenses prior to treatment.

### Cautions

Cautions are important information to prevent possible harm when using Valeda.

- Use only the accessories specified for use with this system. Use of accessories not specified for use with this system may impair performance.
- Only connect a Valeda LumiKey to the system USB port. Connection to other equipment could result in previously unidentified risks to patients, operators or third parties. The manufacturer has implemented measures to control risks related to USB port connections. Changes to the USB port connection can introduce new risks.
- In the event of an emergency in which power to the system must be shut off, either turn off the power switch or unplug the system from the electrical source.
- Do not attempt to disassemble Valeda.

### Precautions

Precautions are important information for alerting the operator to exercise special care necessary for the safe and effective use of Valeda or a measure taken in advance to avert possible harm or misfortune.

- This device has not been tested for use on pregnant women.
- Do not spray or pour cleaning agents directly on Valeda.
- Do not use abrasive agents to clean Valeda.

### Information

The following useful information is related to the operation of Valeda.

- There are light emitting diodes (LEDs) in Valeda emitting at 590, 660 and 850 nm from the light aperture. Valeda has been designed to limit LED output to be within acceptable industry standards (IEC 62471) for eye safety.

### Clinical Data

Information related to the clinical safety and efficacy of Valeda can be found in the LumiThera Report of Prior Investigations and is available upon request.

## 8. Installation

### Overview

Valeda is intended for use in a controlled room environment at ambient temperature and humidity. See Technical Specifications on page 23.

Before unpacking Valeda, inspect the shipping box for damage. If damage is present, notify the freight carrier and ask for an agent to be present while the system is unpacked. Retain all packing materials in their original condition and contact LumiThera.

### Box Contents

The shipping box contains the following:

- Valeda Light Delivery System
- Valeda LumiKey
- Package of Valeda Forehead Rest Protective Barriers
- Package of Valeda Chin Rest Protective Barriers
- Valeda Dust Cover
- Valeda Power Cord(s)
- Valeda Light Delivery System User Manual
- Quick Reference Guide

### Valeda Setup

- Valeda should be used with chairs for the patient and operator and a height-adjustable table to accommodate variations in patient height.
- Remove the system from the shipping box and place on the adjustable table.
- Position Valeda so that the operator can access both the patient and operator sides of the system. The patient must be able to be comfortably seated at the patient side of the system.
- Ensure that Valeda is on a flat, hard surface and that the side vent slots are not blocked.
- Multiple power cords may be in the box. Select the appropriate power cord for your location. Attach the power cord to the unit and connect to a grounded outlet.
- Ensure the power cord is not a tripping hazard.
- Following hot or cold storage (below 15 or above 30°C), allow one hour for the system to reach ambient temperature and humidity before using.

## 9. Operating Instructions

### Operation of Valeda consists of the following steps:

- System setup
- Patient preparation
- Beam alignment procedure
- Treatment procedure
- Cleaning and storage

### Valeda Setup

- Remove the dust cover.
- Verify the system has no obvious damage.
- Confirm the system is plugged into an appropriate electrical outlet and all connections are secure. See Technical Specifications on page 23.
- Turn on the system by moving the power switch to the on position.
- Wait for the system to power on and run through the self-test. When the self-test is complete, the home screen will appear on the touchscreen.

FIGURE 4

To initiate a treatment for a patient, tap the large “+” button (Figure 4).

The next screen shows the instruction for applying the disposable patient protective barriers, one to the forehead rest and one to the chin rest (Figure 5).

Apply the disposable protective barriers to the chin and forehead rest and tap the “Continue” button. →


FIGURE 5

The next screen shows the number of Valeda Treatment Credits remaining. Touch the “Continue” → button (Figure 6).

Refer to the Treatment Credit Section if credits are needed.

FIGURE 6

Select whether one eye or both eyes will be treated (Figure 7).


FIGURE 7

If one eye was selected for treatment, select on the touchscreen which eye is the one to be treated (Figure 8). If both eyes were selected for treatment, this screen is not presented to the operator. The right (OD) eye will always be treated first, followed by the left (OS).


FIGURE 8

The camera will move to the chosen treatment eye and will turn on to allow viewing and positioning of the beam on the patient's eye.

#### **Patient Preparation**

1. Verify that the patient has removed glasses and contact lenses.
2. Position the patient's forehead against the forehead rest. If the patient's forehead is not detected, the system will display "No patient detected."
3. Adjust the chin rest elevation to position the patient's eyes vertically in the centre of the light aperture (Figure 9).
4. Adjust the patient's chin position on the chin rest as shown below (Figure 9).


FIGURE 9

5. Confirm the patient's hair is outside the light aperture.
6. Monitor the patient throughout treatment to ensure proper patient alignment with the system is maintained.

### Beam Alignment Procedure

While watching the camera view on the touchscreen, use the beam alignment joystick on the operator side of Valeda to move the targeting beam (yellow circle) left/right and up/down until the beam fills the eye socket. The picture on the touch screen illustrates the correct position of the beam and the patient's eye to be viewed (Figure 10).

FIGURE 10

Once the beam is aligned, instruct the patient to remain in position until instructed to move.

Press the "Confirmation"  button on the touchscreen (Figure 10).

### Treatment Procedure

The treatment will consist of 4 phases:

1. The treatment initiates with 35 seconds of pulsed yellow and near-infrared (NIR) wavelengths, eyes open.
2. This is followed by 90 seconds of continuous red light wavelength, eyes closed.
3. The treatment will repeat with another 35 seconds of pulsed yellow and NIR wavelengths, eyes open.
4. The treatment will finish with another 90 seconds of continuous red light wavelength, eyes closed.

The operator will instruct the patient to keep their eyes open during phases 1 and 3 and then to keep their eyes closed for phases 2 and 4.

The touchscreen will show the Treatment Ready screen.


FIGURE 11

Verify that all the information on this screen is correct, then touch the “Confirmation”  button (Figure 11).

The patient is instructed to keep the eye open for phase 1 (Figure 12).


FIGURE 12

Press the **green** activation button on the system to initiate treatment. A progress screen will be shown. A countdown timer is displayed while the progress bar fills from left to right (Figure 13).


FIGURE 13

Valeda will pause after completion of each phase. For example, after the first wavelength treatment phase is complete, the screen will provide instructions for the patient's eyelid position (open or closed) and display a live video for confirmation of the eyelid position (Figure 14).


FIGURE 14

The patient is instructed to keep their eyes closed for phase 2.

To proceed with the next phase, press the **green** activation button on the system. A countdown timer is displayed while the progress bar fills from left to right (Figure 15).


FIGURE 15

The next steps to complete phases 3 and 4 are identical to the steps detailed for phases 1 and 2.

If the patient's second eye is to be treated, the system will display the following screen after the 4 phases of treatment for the first eye are completed (Figure 16).


FIGURE 16

Touching the "Confirmation"  button will return to the Beam Alignment Procedure for treatment of the second eye.

Upon completion of the treatment session, the finished screen will be displayed (Figure 17).

Instruct the patient to pull away from the system.

Remove and dispose of the forehead and chin rest protective barriers.

FIGURE 17

To continue patient treatments, touch the "Home"  button.

If all patients have been treated for the day, switch off the power to Valeda.

### Treatment Pause

A treatment phase may be paused at any time by pressing the **green** activation button on the system, or by the patient pulling back from the system (Figure 18).


FIGURE 18

To resume treatment, press the **green** activation button on the system.

Check the alignment and state of the eyelid on the video screen (Figure 19).


FIGURE 19

Press the **green** activation button to restart the treatment at the time it was paused.

### Treatment Cancellation

While the treatment is paused, the treatment can be cancelled by pressing the **red** button on the system; the message below will be displayed (Figure 20). If a treatment is cancelled, the patient has not completed a full treatment. The physician and patient can determine whether the treatment may need to be rescheduled.

FIGURE 20

### Clean and Store System

Remove and dispose of the disposable protective barriers on the forehead rest and chin rest. If needed, clean the system per cleaning instructions provided in this User Manual (page 20). Confirm the system is off and place the dust cover over the system.

## 10. Valeda Treatment Credits

### Valeda Will Not Operate Without Treatment Credits

- Treatment credits are purchased on the LumiThera e-commerce store at [www.lumithera.com](http://www.lumithera.com).
- LumiThera will provide an electronic LumiFile.
- Copy the provided electronic file to your LumiKey.
- Insert the LumiKey into the USB port on the front of Valeda.
- Valeda will automatically read the LumiKey and add the appropriate number of credits to your system.

### Grace Period Treatment Credits

Grace credits are to be used in the event that online purchasing is unsuccessful. Contact LumiThera customer care to purchase grace credits.

To enable grace period treatment credits:

- From the main menu, press the “Information”  button.
- At the top right corner of the information screen, press the “+” button.
- LumiThera will provide an alphanumeric code that is specific to your system. Grace period codes must be loaded onto the system within 24 hours.
- Enter the code into Valeda and press the “Continue”  button to proceed (Figure 21).


FIGURE 21

If the grace period code was accepted, a confirmation screen will be displayed with the total credits available (Figure 22).


FIGURE 22

If the grace period code was rejected, the following screen will appear (Figure 23).


FIGURE 23

## 11. Valeda Settings

Press the “Home”  button. The Settings Menu may be reached by pressing the “Settings”  button on the display (Figure 24).

FIGURE 24

Selecting the “Language” button provides a listing of languages available on your system (Figure 25). Not all systems will include multiple language options.

FIGURE 25

Pressing the “Information”  button will display key system information including treatment credits remaining, maintenance dates and software versions (Figure 26).

FIGURE 26

## 12. Cleaning

- The patient light aperture may be cleaned with a lens cleaning cloth or optical glass cleaner.
- Patient contact areas may be cleaned with isopropyl alcohol wipes.
- The exterior of Valeda, not including the patient light aperture, may be cleaned with a mild soap solution.
- Do not spray or spill any fluid onto the unit. Use a damp or dry cloth for cleaning.
- Contact LumiThera if the light aperture appears contaminated.

## 13. Disposal

### Protective Barrier Disposal

Follow your practice’s procedure, industry guidelines or government regulations applicable to the disposal of protective barriers.

### System

Disposal of electronic equipment may be regulated in your community. Follow local regulations for disposal of electronic equipment.

## 14. Maintenance

- Each Valeda has been manufactured to specification.
- Operator maintenance is not required.
- All service and maintenance will be provided by LumiThera.
- If there is a problem, or it is suspected that Valeda requires maintenance, contact LumiThera.

**CAUTION** – Do not attempt to disassemble the system.

There are no operator serviceable or repairable components. Do not attempt to take the assembly apart.

## 15. Troubleshooting

In the event a system fault is detected by Valeda, a fault message will be displayed on the screen. If possible, the display screen will identify the fault (Figure 27).


FIGURE 27

For each system fault, perform the following:

- Record the fault message.
- Rebooting the system may allow a treatment to proceed.
- Contact LumiThera for assistance.
- Report all system issues to LumiThera.

## Troubleshooting Guide for Errors and Alarms

Display Message	Description	Additional Instruction
Overheated LT-166	The LEDs became too hot during operation.	This can be caused by the ambient temperature in the room being too high. The system must not be operated in a room where the ambient temperature is above 30°C. The system may overheat if the cooling fan vents are blocked. Make sure the system is on a flat, hard surface and nothing is blocking the vent slots.
Patient sensor failed LT-141	The patient sensor is mounted in the system forehead rest.	The system must be repaired by a LumiThera service technician.
Eye select mirror stalled LT-133	The internal motor that positions the light beam for the correct eye has become stuck.	The system must be repaired by a LumiThera service technician.
LED failed LT-117	One of the LEDs is not working.	The system must be repaired by a LumiThera service technician.
Light engine communication error LT-190	The touchscreen computer cannot communicate with the light engine.	The system must be repaired by a LumiThera service technician.
Cooling fan error LT-232	The fan in the system has a sensor that detects if the fan is not rotating.	The system must be repaired by a LumiThera service technician.
Camera error LT-208	The camera used for targeting has failed.	The system must be repaired by a LumiThera service technician.
Light engine failed LT-109	The light engine has built-in safety checks. If any of these safety checks fail, this message is displayed and the system cannot be operated.	The system must be repaired by a LumiThera service technician.
Start or stop button failed LT-158	One of the buttons on the front of the system is not working correctly.	The system must be repaired by a LumiThera service technician.

## 16. Technical Specifications

Parameter	Specification
Size	530 mm height x 300 mm width x 330 mm depth
Weight	10.8 kg
Electrical Power Input	100 – 240 VAC, 50/60 Hz, 1.0 – 0.5 A
Fuse Rating	Voltage: 250 V Current: 2 A Operating speed: Slow Size: 5 x 20 mm IEC Fuse Code: T2AL250V
Light Sources	Light emitting diodes (LEDs)
Light Emission	590 nm output: 5 mW/cm <sup>2</sup> 660 nm output: 65 mW/cm <sup>2</sup> 850 nm output: 8 mW/cm <sup>2</sup>
Beam Diameter	30 mm (nominal) at treatment plane
Treatment Exposure Time	A total of 250 seconds (4 minutes 10 seconds). There are 4 phases: Phase 1 – 35 seconds, patient's eyes open Phase 2 – 90 seconds, patient's eyes closed Phase 3 – 35 seconds, patient's eyes open Phase 4 – 90 seconds, patient's eyes closed
Operating Environment	Temperature: 15 – 30°C Humidity: 15 – 90% (non-condensing) Altitude: Up to 3,000 m
On-site Storage Environment	Temperature: 15 – 30°C Humidity: 15 – 90% (non-condensing) Altitude: Up to 3,000 m
Transportation Environment	Temperature: -18 - 60°C Humidity: 15 - 90% (non-condensing) Pressure: 700 - 1060 hPa
Ingress Rating	EC 60529: IP2X
Electrical Safety Classification	IEC 60601-1 Type BF applied part
Electromagnetic Compatibility	IEC 60601-1-2
Photobiological Safety	IEC 62471 Risk Group: Exempt

## Technical Specifications Continued

Parameter	Specification
Valeda LumiKey	USB 2.0-compliant flash drive Compatible with Microsoft Windows USB Mass Storage class drivers Can be used with Apple Mac computers that have USB A ports

## 17. Glossary

Acronym/Term	Definition
AC	Alternating current
A/m	Amperes per metre
Aperture	An opening through which light energy is emitted
ESD	Electrostatic discharge
GHz	Gigahertz
Hz	Hertz
IEC	International Electrotechnical Commission
Irradiance	Radiant power per unit area (radiant power divided by the area of the beam); i.e., the quotient of the radiant flux ( $d\Phi$ ) incident on an element of a surface by the area $dA$ of that element (SI units = $W/m^2$ )
kg	Kilograms
kV	Kilovolts
kHz	Kilohertz
LED Power	The time rate of flow of radiant energy, measured in joules per second or watts
LED	Light emitting diode
MHz	Megahertz
m	Metre
mm	Millimetre, $1 \times 10^{-3}$ of a metre
Near Infrared Radiation	Light that is of slightly longer wavelength than the human eye can see

## Glossary Continued

Acronym/Term	Definition
nm	Nanometre, $1 \times 10^{-9}$ of a metre
OD	Oculus dexter, the patient's right eye
OS	Oculus sinister, the patient's left eye
P	As defined in section 18
PBM	Photobiomodulation
Pulsed	LED output rapidly transitioning between on and off states
RF	Radio frequency
$\sqrt{P}$	Square root of P (as defined in section 18)
V/m	Volts per metre
VAC	Volts alternating current
VRMS	Volts root mean square
W	Watts

## 18. Guidance and Manufacturer's Declarations

### Electromagnetic Emissions for All Equipment and Systems

Valeda is intended for use in the electromagnetic environment specified below. The operator of Valeda should ensure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF Emissions CISPR 11	Class A	The system is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Power Line Fluctuations (Flicker) EN 61000-3-3	Complies	

### Electromagnetic Immunity for All Equipment and Systems

Valeda is intended for use in the electromagnetic environment specified below. The operator of Valeda should ensure that it is used in such an environment.

Immunity Test	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 15 kV air ± 8 kV contact	Floors should be made of wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical Fast Transient/Burst IEC 61000-4-4	2 kV at AC power lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge Transient IEC 61000-4-5	2 kV line-to-earth at AC power ports  1 kV line-to-line at AC power ports	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips and Interruptions IEC 61000-4-11	- 0%, 0.5 Cycle - 0%, 1 Cycle - 70%, 25 Cycles - 0%, 250 Cycles (interruption)	Mains power quality should be that of a typical commercial or hospital environment.
Power Frequency Magnetic Field IEC 61000-4-8	30 A/m; 50/60 Hz	Power frequency magnetic fields should be that of a typical commercial or hospital environment.

### Electromagnetic Immunity for Equipment and Systems That Are Not Life-Supporting

Valeda is intended for use in the electromagnetic environment specified below. The operator of Valeda should ensure that it is used in such an environment.

Immunity Test	Compliance Level	Electromagnetic Environment – Guidance
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	<p>Portable and mobile RF communications equipment should be used no closer to any part of the system than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p><b>Recommended separation distance</b>  <math>d = (1.2)(\sqrt{P})</math></p> <p><math>d = (1.2)(\sqrt{P})</math> 80 to 800 MHz</p> <p><math>d = (2.3)(\sqrt{P})</math> 800 MHz to 2.5 GHz</p> <p>Where P is the maximum power output rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.</p> <p>Interference may occur in the vicinity of equipment containing a transmitter.</p>

### Recommended Separation Distances for Equipment and Systems That Are Not Life-Supporting

Valeda is intended for use in an electromagnetic environment in which radiated disturbances are controlled. The operator of Valeda can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment and the system as recommended below, according to the maximum output power of the communications equipment.

Maximum Output Power (watts)	Separation (metres) 150 kHz to 80 MHz $d = (1.2)(\sqrt{P})$	Separation (metres) 80 MHz to 800 MHz $d = (1.2)(\sqrt{P})$	Separation (metres) 800 MHz to 2.5 GHz $d = (2.3)(\sqrt{P})$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

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LumiThera, Inc.