

THIS IS WHAT THE FUTURE LOOKS LIKE

Castroviejo

Corneoscleral Scissors small blades K4-1000 right (shown) K4-1010 left



Castroviejo

Corneoscleral Scissors medium blades K4-1020 right (shown) K4-1030 left



Troutman-Castroviejo

Corneoscleral Scissors with stop, small blades K4-1200 right (shown) K4-1210 left







Troutman-Castroviejo

Corneoscleral Scissors with stop, medium blades K4-1220 right K4-1230 left (shown)



Troutman-Katzin

Corneal Transplant Scissors with stop strong curve, small blades K4-2000 right K4-2010 left (shown)



Rich

Corneal Transplant Scissors extremely delicate, strongly curved blades For trimming the recipient corneal bed and donor button. Can also be used for open sky anterior capsulotomy. K4-2014 right K4-2015 left (shown)





ΟΡΗΤΑΙΜΟ

Terry

Micro Corneal Scissors round handle for meticulous control K4-2016 right (shown) K4-2017 left



Holland

DALK Scissors spatulated tip K4-2018 right (shown) K4-2019 left



Katzin

Corneal Transplant Scissors strong curve, medium blades K4-2020 right K4-2030 left (shown)









Castroviejo

Universal Corneal Scissors curved, blunt tips K4-2200 small blades K4-2220 medium blades K4-2240 large blades K4-2260 strongly curved large blades



Castroviejo

Keratoplasty Scissors angled to side blunt tips, medium blades K4-2850



McPherson-Westcott

Conjunctival Scissors curved blunt tips, small blades K4-3000









Sarkisian-Westcott

Scissors small blades with nine 1mm landmarks laser etched on one blade For intra-operative cutting and measuring with a single instrument. K4-3002



Westcott

Tenotomy Scissors blunt tips, medium blades curved K4-3004

Westcott

Tenotomy Scissors blunt tips, standard blades curved K4-3100 right (shown) K4-3110 left







Tibolt

Punctal Scissors one conically shaped blade with micro grooves and textured surface For simple "Two-snip Punctoplasty". K4-3130



Mini Westcott

Tenotomy Scissors 4mm long curved blades blunt tips K4-3200



McPherson-Westcott

Stitch Scissors very sharp pointed tips small blades curved K4-4000







Westcott

Stitch Scissors very sharp pointed tips standard blades curved K4-4100



Jaffe

Stitch Scissors very sharp pointed tips medium blades curved K4-4150



Katena

Micro Tying Scissors For tying and cutting 9-0 to 11-0 sutures. 2.75mm long tying platform Not to be used as a needle driver. K4-4200









Iris Scissors sharp tips, small blades K4-4400 straight K4-4410 curved



Fine Stitch Scissors sharp tips, medium blades K4-4500 straight K4-4510 curved



Barraquer

Iris Scissors blunt tips, 7mm long blades K4-4900



ΟΡΗΤΑҍΜΟ



Capsulotomy Scissors

Capsulotomy Scissors


Vannas 5mm long blades, sharp tips K4-5000 straight K4-5010 curved K4-5020 angled forward K4-5030 angled to side

Gills-Vannas

7mm long blades, sharp tips K4-5050 straight K4-5060 curved K4-5070 angled forward

Holland DALK Scissors 4.5mm angled blades K4-5075 blunted tips

McPherson-Vannas 8mm long blades, sharp tips K4-5100 curved



K4-5075

K4-5070


Capsulotomy Scissors

Katena-Vannas 7mm extra thin long blades, sharp tips K4-5110 Katena-Vannas, straight K4-5111 Katena-Vannas, curved K4-5112 Gills-Vannas, angled forward

10mm extra thin long blades, sharp tips K4-5120 Katena-Vannas, straight K4-5121 Katena-Vannas, curved K4-5122 Stern-Gills, angled forward K4-5123 Uribe-Stern, curved,

10mm long angled blades, curved tipsK4-5125Gills, sharp pointsK4-5126Uthoff-Gills, blunt points

angled, forward





K4-5060



K4-5100

Vitrectomy

Scissors 20 gauge shank Stainless steel construction with serrated plastic handle. Autoclavable.



Scissors horizontal opening angled 30° K4-5210



Scissors vertical opening angled 60° K4-5240



Scissors vertical opening angled 80° K4-5250



Katena

"Squeeze Handle" Scissors These "squeeze handle" instruments feature titanium handles to reduce overall weight, and front ends of hardened stainless steel for secure cutting.



Scissors horizontal, straight K4-5300 20ga.

Scissors horizontal, curved **K4-5310** 20ga. **K4-5423** 23ga. **K4-5425** 25ga.



Scissors vertical, 40° K4-5330 20ga. Scissors vertical K4-5433 23ga.



Scissors horizontal, angled 45° K4-5320 20ga.



Scissors vertical, 80° K4-5340 20ga.









Alio

MICS Scissors 23 gauge curved shaft with horizontal micro blades For cutting delicate membranes. K4-5351



Scimeca

Lens Loop Amputator For cutting polypropylene or PMMA lens haptics in the anterior or posterior chamber. K4-5550 with guard



OPHTAĿMO

IOL Cutters

Donnenfeld

Micro IOL Cutter 4mm long blades micro notches in blades for traction For cutting soft IOLs through a 1.5mm incision. K4-5560





Stolte

IOL Cutter titanium slotted cover to retain soft IOL during cutting process For cutting soft IOLs. K4-5565





Katena

IOL Cutter 6mm long serrated blades For cutting soft IOLs. K4-5571



Eye Scissors	
pointed tips	
small	
K4-7200	straight
K4-7210	curved

Eye Scissors

pointed tips standard K4-7400 straight K4-7410 curved

Eye Scissors

large rings pointed tips standard K4-7440 straight K4-7450 curved



10

0

DESIGNED FOR SIGHT

Stitch Scissors ribbon type needle points curved K4-8000




Stitch Scissors needle points light curve K4-8100











ΟΡΗΤΑμΜΟ

Strabismus Scissors

ribbon type blunt tips K4-9020 straight K4-9030 curved



Metzenbaum

Dissecting Scissors gently curved blades blunt tips 5¾" (14.5cm) long K4-9040

















O SURGICAL

• CONSUMABLES

• MEDICAL DEVICE

OPHTALMO | Industriepark De Bruwaan 5 | B-9700 Oudenaarde | T + 32(0) 5531 3038 | ophtalmoservice.be